



Cleaning Tool


For automated cleaning




For manual cleaning



For drying



Cleaning Tube



Primado2
Slim Attachment 300



NAKANISHI INC.
700 Shimohinata, Kanuma Tochigi 322-8666, Japan
TEL : +81 289 64 3380 FAX : +81 289 62 5636

Slim Attachment 300 80,000min⁻¹

Slim Attachment 300 is designed with new features to deliver improved visibility, smooth operation, safety, and stable rotation performance.

Variable Exposure

- Adjustable exposure lengths of 7/9/11/13/15/17mm

Non-Slip Surface

- Improved grip with surface etching

Quick Joint

- Easy attachment/detachment with Quick Joint
- 16 different handpiece combinations are possible with the straight/angled hubs and variety of tubes

Smooth Connection


- Quick Joint allows for an easy attachment-motor connection

Accessory

- Two types of accessories can slide onto the attachment

Tapered Tip

- Ensures better visibility and access to surgical anatomy



Compatible with Primado 2 Control Unit and Slim Motor



order code
P200-CU-120
Primado2 Control Unit for 120V



order code
P200-CU-230
Primado2 Control Unit for 230V



order code
P200-SMH-HS
Primado2 Slim Motor/HS
· Hand Control Motor
· Speed : Max. 80,000 min⁻¹




order code
FC-73/FC-74
Primado2 Single Foot Control
Primado2 Multi Foot Control



order code
P200-SMH
Primado2 Slim Motor
· Speed : Max. 80,000 min⁻¹


Operability

Straight attachment added to the product lineup.



Reduce hand fatigue with the Straight Hub.



By using a Straight Hub, surgeons can easily access a narrow space by simple rotation.



Non-slip surface etching improves grip.

Accessibility

Better accessibility with the 6 exposure lengths adjusting mechanism. The length can be adjusted with the Accessory attached.




Visibility

Visibility is improved with the tapered tip and variable exposure lengths.

exposure length 11mm



exposure length 17mm



Lineup

Slim Attachment Hub



order code
P300-1AHS
Slim Attachment Hub/Straight



order code
P300-1AHA
Slim Attachment Hub/Angled

Slim Tube



order code
P300-1T110
Slim Tube 110



order code
P300-1T110-C
Slim Tube 110/C



order code
P300-1T110-PB
Slim Tube 110/PB



order code
P300-1T130
Slim Tube 130



order code
P300-1T130-C
Slim Tube 130/C



order code
P300-1T130-PB
Slim Tube 130/PB



order code
P300-1T170
Slim Tube 170



order code
P300-1T240
Slim Tube 240

Accessory



order code
P300-1T-B
Slim Tube Beak



order code
P300-1T-H
Slim Tube Hood

*Cannot be used with the Irrigation Nozzle.

Irrigation Nozzle



order code
P200-IN-1A155
For P300-1T110/P300-1T110-C/P300-1T110-PB

order code
P200-IN-1A175
For P300-1T130/P300-1T130-C/P300-1T130-PB



Cutting Accessories for Slim Attachment 300

Cutting Accessories							Head size ● (mm) Actual size	Accessory	Slim Tube	Slim Attachment Hub
R Fluted Burs	R Diamond Burs	R Coarse Diamond Burs	MH Diamond Burs	MH Coarse Diamond Burs	RT Diamond Burs	RT Coarse Diamond Burs				
order code	order code	order code	order code	order code	order code	order code				
—	PDS-1RD110-20	PDS-1RCD110-20	PDS-1MHD110-20	—	—	—	● 2.0	 order code P300-1T-B Slim Tube Beak	110 mm	 order code P300-1T110 Slim Tube 110  order code P300-1T110-C Slim Tube 110/C  order code P300-1T110-PB Slim Tube 110/PB  order code P300-1AHS Slim Attachment Hub/ Straight
PDS-1RF110-30	PDS-1RD110-30	PDS-1RCD110-30	PDS-1MHD110-30	PDS-1MHCD110-30	—	—	● 3.0	 order code P300-1T-H Slim Tube Hood		
PDS-1RF110-40	PDS-1RD110-40	PDS-1RCD110-40	PDS-1MHD110-40	PDS-1MHCD110-40	—	—	● 4.0			
—	—	—	—	—	PDS-1RTD110-45	PDS-1RTCD110-45	● 4.5	—		
—	—	PDS-1RCD110-50	—	—	—	—	● 5.0			
—	PDS-1RD130-20	PDS-1RCD130-20	PDS-1MHD130-20	—	—	—	● 2.0	 order code P300-1T-B Slim Tube Beak	130 mm	
PDS-1RF130-30	PDS-1RD130-30	PDS-1RCD130-30	PDS-1MHD130-30	PDS-1MHCD130-30	—	—	● 3.0	 order code P300-1T-H Slim Tube Hood	 order code P300-1T130 Slim Tube 130	
PDS-1RF130-40	PDS-1RD130-40	PDS-1RCD130-40	PDS-1MHD130-40	PDS-1MHCD130-40	—	—	● 4.0		 order code P300-1T130-C Slim Tube 130/C	
—	—	—	—	—	PDS-1RTD130-45	PDS-1RTCD130-45	● 4.5	—	 order code P300-1T130-PB Slim Tube 130/PB	
—	—	PDS-1RCD130-50	—	—	—	—	● 5.0			
PDS-1RF170-35	PDS-1RD170-35	PDS-1RCD170-35	—	—	—	—	● 3.5		170 mm	
PDS-1RF170-40	PDS-1RD170-40	PDS-1RCD170-40	—	—	—	—	● 4.0		 order code P300-1T170 Slim Tube 170	
—	PDS-1RD170-50	PDS-1RCD170-50	—	—	—	—	● 5.0			
—	PDS-1RD240-30	PDS-1RCD240-30	—	—	—	—	● 3.0	 order code P300-1T-B Slim Tube Beak	240mm	
PDS-1RF240-40	PDS-1RD240-40	PDS-1RCD240-40	—	—	—	—	● 4.0	 order code P300-1T-H Slim Tube Hood		
—	PDS-1RD240-45	PDS-1RCD240-45	—	—	—	—	● 4.5		 order code P300-1T240 Slim Tube 240	
—	PDS-1RD240-50	PDS-1RCD240-50	—	—	—	—	● 5.0			